Evaluation of Corrosion Effects on Surgical Instruments



conference of AFS

Dr. Gerhard Kirmse, Tuttlingen, Germany

Conclusions

- Corrosion represent a common issue in hospitals
 Different evaluations of risk situation triggered further research
- · high refurbishment costs can be avoided by process optimization


RISKS

- · in the literature some case studies indicate a relation between post operative
- wound healing disorders and instrument residue randomly selected corroded instrument analysis shows significant amounts
- amounts of protein range to critical amounts close to 100µg


ROOT CAUSES

· Water Quality can be identified in standardized simulations as a root cause for corrosion => demineralized water according to EN285 is a prerequisite to avoid corrosion

• other identified root causes are lack of lubrication (friction corrosion) and cleaning performance



Laboratory Analysis of Surface Changes



Root Causes of Surface Changes (most frequent ,out of over 100 process optimizations)



Contact: Dr. Gerhard Kirmse (gerhard.kirmse@aesculap.de) Dr. Kirmse is an employee of Aesculap AG, Tuttlingen and member of the AKI (Working Group Reprocessing)

Literature Analysis of Instruments as cause of Infections

- Limited literature situation
- · Research typically done upon raise of post-operative wound infections in some cases links to instrument residue/ corrosion were found
- · however no full proof of root causes cause is possible

Examples

Seattle Children's warns of potential infection risk due to "debris" on instruments

 Dancer (Scotland) reported about conatminated instruments after faulty packaging and sterilization ektionen im Zusam


nenhang mit kontaminierten chirurgischen Instru hanie J. Dani 3:2: 134-135 Infections with Pseudomonas Aeroginosa were reported and linked to improperly

cleaned shaver handpieces in the US in 2012 The hid ns. Iwatch News, 22.02.2012

Mamalis described metal oxide as one root cause of the Toxic-Anterior Segment

Syndroms (TASS) in eye surgery,


Examples



Results:

. In 24 of 34 analyzed instruments, the elementary carbon content of residue was more than 10% • Co-Elements like Nitrogen, Sulfur, Potassium frequently point

towards organic / patient residue



further information. http://www.a-k-i.org "Proper Maintenance of Instruments" http://www.insightsforexcellence-aesculap.com